

Medical Gas Area Alarm Medipoint 26 Digital (MP26D) Local Alarm

SPECIFICATION

Medical Gas Area Alarm

Each medical gas area alarm panel shall be capable of monitoring 6 medical gas services by means of pressure sensors, which detect deviations from the normal operating limits of either pressure or medical vacuum. The medical gas area alarm shall fully comply with the requirements of HTM2022, HTM 02-01 and BS EN ISO 7396-1.

The cover, backbox and bezel (if required) shall be polyester powder coated in a RAL9010 30% gloss finish. A single tamperproof fastener shall be used to gain access to the hinged door. The hinge shall operate through a minimum of 120° to provide adequate access.


Antimicrobial Additive

User accessible parts such as the screen protector, key pad fascia shall include a silver antimicrobial additive for inherent antimicrobial protection.


Alarm Status Screen

Each gas service shall be displayed by coloured 7" LCD touch screen to display the pipeline pressure and status for normal and fault conditions in line with HTM02-01 and HTM2022 requirements. Screen graphics are available in normal and dark mode settings for reduced brightness.

Typical normal mode view



Typical dark mode view



Gas Service Status Badges

The gas status bar shall indicate green steady for normal condition, and change to red flashing for fault conditions. Loss of communication with the gas sensor is indicated by changing the pressure reading to "ERR" and gas ID bar to red steady. The gas pressure is displayed for each service along with a vertical bar showing the alarm set points relative to the live pressure reading.

Normal	Pressure Fault	Sensor Fault
OXYGEN/NITROUS OXIDE 4.2 Bar NORMAL	OXYGEN/NITROUS OXIDE 3.4 Bar LOW PRESSURE	OXYGEN/NITROUS OXIDE 5.2 Bar HIGH PRESSURE
		OXYGEN/NITROUS OXIDE ERR FAULT

Alternative Sensor Configuration

Each input can be configured for many common 4-20 mA sensors inputs including alarm trigger points from the on screen setup menus.

Typical examples of customised sensors

MEDICAL O2 FLOW	PRESSURE DEW POINT	CO CONCENTRATION
154 l/min NORMAL	-49 °C NORMAL	22 PPM NORMAL

Audible Alarm and Mute function

An audible warning shall sound simultaneously with any failure indication and a mute facility shall be provided. Following a mute selection the audible will resound after approximately 15 minutes, or shall operate simultaneously should a further alarm condition occur. A "Mute" switch shall be provided inside the panel for use during any maintenance resulting in prolonged pipeline or plant shutdown. This facility shall automatically reset when the gas service returns to normal.

Icon	Description
Mute Button	 Press to activate to silence the audible alarm.

Test Function

The alarm panel shall have a 'Test' facility to prove the integrity of the internal circuits, gas status and audible warning.

Icon	Description
Test Button	 Press to activate the alarm test button.



Editable Alarm Service Location Legend

The top strip on the screen display has a field for displaying the details for the department location served by the gases that can be edited through the setup menus.

Main Screen Interactive Icons


There are interactable icons for accessing the alarm information screen, alarm logs and setup menus. Alarm logs and setup menus are password protected to prevent unauthorised access to these screens.

Icon	Description
	Info Icon. Takes you to the info screen
	Logs Icon. Takes you to the logs screens.
	Settings Icon. Takes you to the settings menu screens.

Interactive Setup Menus

From accessing the on screen setup menus the alarm can be configured to match the hospitals medical gas alarm system specifications.


Typical menu setup screen



The screenshot shows a 'LOCATION' setup screen. On the left, there's a display for 'OXYGEN' showing '4.4 Bar' and a 'NORMAL' status bar. Below this is a 'ZERO' button. On the right, there's a 'SELECT GAS' button, a 'GAS ID' dropdown set to '00', an 'INPUT' dropdown set to '1', and two numeric input fields for 'LOW T/H' (3.6) and 'HIGH T/H' (4.9). At the bottom right are 'X' and 'checkmark' buttons.

On screen key pads allow for simple editing of text fields. Including main keyboard, alternative keyboard and numpad.

Typical on screen keyboard



The screenshot shows a virtual keyboard with blue buttons. The layout includes a top row (Q, W, E, R, T, Y, U, I, O, P), a second row (A, S, D, F, G, H, J, K, L), a third row (Z, X, C, V, B, N, M, left arrow), and a bottom row with an 'ALT' button, a long spacebar, and an 'OK' button.

Power on and System Fault LED Indicators

Each alarm shall provide a green LED to indicate that electrical power is available at the panel and a red LED to indicate 'System Alarm'. In the event of a system fault or electrical power supply failure the 'System Alarm' LED shall illuminate (flashing).

Icon	Description
Power LED	 Green, illuminated - Power On
	 Black, de-energised - Power Off
System Fault	 Black, de-energised - System Normal
	 Red, illuminated - System Fault

Sensor or Communication Fault Monitoring

Line contact monitoring circuits shall be provided to constantly monitor the integrity of the input sensors and interconnecting wiring. In the event of any fault the line contact monitoring circuits shall initiate the specific gas service failure indication, a 'System Alarm' indication and an audible warning. Further aids to fault diagnosis shall be provided by means of varying flashing rates whilst operating the 'Test' switch.

Alarm System Networking.

A simple data connection shall be provided to allow connection, enabling the visual and audible alarm signals to be repeated at other locations within a department.

Pressure Sensors

Pressure sensor scale shall include negative and positive pressure scale, allowing one sensor to cover all gas and vacuum supply monitoring. Electrical connectors shall be designed for frequent disassembly. Spade connectors are not acceptable. Pressure switches shall include a 1/4" BSP threaded pipeline connection.



Relay Output

The alarm includes a relay output connection to allow the overall status of the alarm to other alarms or building management systems. Any fault condition on the alarm will change the status of the relay output.

Remote Audible

The alarm includes a connection for a remote buzzer. Any condition that activates the alarm's internal audible warning will also activate the remote buzzer. Mute functions will silence the remote buzzer as per the internal audible warning.

Data Logging

Data logs can be accessed through password protected screens to review up to 1000 logged events. Data logging includes a main screen from all events, plus filtered screens to see each gas service separately.

Typical data logging screen

LOCATION									
Log Viewer: Service 1 Gas O2 (Showing Entries 1 to 15 of 35)									
ID	Activated	Gas	Ch.	Status	Location	Cleared			
1	12-07-2020 17:54:45	O2	0	High Pressure	WARD 23	12-07-2020 17:55:45			
2	12-07-2020 17:50:11	O2	0	Low Pressure	WARD 23	12-07-2020 17:51:11			
3	12-07-2020 17:47:17	O2	0	High Pressure	WARD 23	12-07-2020 17:48:17			
4	12-07-2020 17:37:15	O2	0	Low Pressure	WARD 23	12-07-2020 17:38:15			
5	12-07-2020 17:37:58	O2	0	High Pressure	WARD 23	12-07-2020 17:38:58			
6	12-07-2020 17:37:56	O2	0	Low Pressure	WARD 23	12-07-2020 17:38:56			
7	12-07-2020 17:36:17	O2	0	High Pressure	WARD 23	12-07-2020 17:37:17			
8	12-07-2020 17:36:03	O2	0	Low Pressure	WARD 23	12-07-2020 17:37:03			
9	12-07-2020 17:36:01	O2	0	High Pressure	WARD 23	12-07-2020 17:37:01			
10	12-07-2020 17:35:34	O2	0	Low Pressure	WARD 23	12-07-2020 17:36:34			
11	12-07-2020 17:34:08	O2	0	High Pressure	WARD 23	12-07-2020 17:35:08			
12	12-07-2020 17:33:58	O2	0	Low Pressure	WARD 23	12-07-2020 17:34:58			
13	12-07-2020 17:33:57	O2	0	High Pressure	WARD 23	12-07-2020 17:34:57			
14	12-07-2020 17:33:55	O2	0	Low Pressure	WARD 23	12-07-2020 17:34:55			
15	12-07-2020 17:33:47	O2	0	High Pressure	WARD 23	12-07-2020 17:34:47			

Power source

Mains operated using 110 V-230 V, 50/60 Hz, alternating current.

Cable size:

Frequency	Voltage	Minimum Cable Size
50 Hz	230 V	3 x 1.5 mm ²
60 Hz	110 V	3 3 AWG14

- Current requirements - 3.0 amps
- Type of protection against electric shock.
- Class 1 (Mains supplied equipment using a protected earth).
- Relative Humidity 90% max.
- Altitude up to 2000 m
- Pollution Degree 2


Mode of operation:

- Indoor use.
- Continuous (equipment may be left switched on indefinitely).

CE Marking


The standard range of BeaconMedæS Medipoint 26 Medical Gas Area Alarms are 'CE' marked with approval from a notified body (more detailed information available on request).

Part Numbers


Part Number	Description	Image
Main Product group		
8102371400	Medipoint 26 Digital Alarm (MP26D)	
4233402315	Medical Gas Alarm 1st Fix backbox	
4233402316	Medical Gas Alarm Bezel	
4233400417	Pressure & Vacuum Sensor	
Related Products		
4233402317	Minimum leak tee adaptor 15 mm	
4233402318	Minimum leak tee adaptor 22 mm	
4233402319	Minimum leak tee adaptor 28 mm	
4233402320	Minimum leak tee adaptor 35 mm	
4233402321	Minimum leak tee adaptor 42 mm	
4233402322	Minimum leak tee adaptor 54 mm	
4233402323	Connection Kit c/w 2mt tube - for remote pressure switch connection	



Installation Details



Wiring Terminal Connections



Ref	Specification	Description
AUD	Remote Audible	Connection for remote buzzer for audible alarm.
SLAVE	N/O Volt free Relay output	Normally open relay output triggered from any alarm fault condition to to another alarm or Building Management System.
G1 to G6	4-20 mA input	Gas sensor or alternative 4-20 mA transmitter input connection
IN	RS-485 input	Medipoint alarm system input communication connection
OUT	RS-485 output	Medipoint alarm system output plus cable screen communication connection
MAINS INPUT	110 V-230 V 50/60 Hz	Electrical power supply connection

